

# Capillary Sequencing (CS)



## Oligonucleotides Synthesis

- Normal Oligos
- Modified Oligos
- Single-labelled Probes
- Double-labelled Probes



## Capillary Sequencing

- Standard Sequencing
- Ready-to-load Sequencing
- PCR Clean-up Sequencing
- Clone Sequencing
- Primer walking Sequencing



## Next-Generation Sequencing

- Genome Sequencing
- Transcriptome Sequencing
- Epigenome Sequencing
- Bioinformatics Analysis



## Customized Project

- Plasmid Preparation
- Nucleic Acid Extraction
- Molecular Cloning
- Bacterial 16sRNA Identification
- Mass Spectrometry

## Select from our wide range of universal primers for your sequencing work

M13 Forward (-20)  
M13-40  
M13 Reverse (-24)  
M13/pUC Forward  
M13/pUC Reverse  
T7  
T7 Promoter  
T7 Minimal Primer (19-mer)  
T7term  
Sp6  
T3  
CMV Forward  
pENTR-F  
pENTR-R  
pGEX 3'  
pGEX 5'  
RVprimer3  
GLprimer2  
AOX1 Reverse  
AOX1 Forward  
BGH-rev  
pET Upstream primer  
GFP-F  
GFP-R  
pBABE 3'  
pBABE 5'  
pTrcHis Forward  
pTrcHis Reverse

5' - GTA AAA CGA CGG CCA GT -3'  
5' - GTT TTC CCA GTC ACG AC - 3'  
5' - GGA AAC AGC TAT GAC CAT G - 3'  
5' - CCC AGT CAC GAC GTT GTA AAA CG -3'  
5' - AGC GGA TAA CAA TTT CAC ACA GG -3'  
5' - GTA ATA CGA CTC ACT ATA GGG C - 3'  
5' - TAA TAC GAC TCA CTA TAG GG - 3'  
5' - TAA TAC GAC TCA CTA TAG G - 3'  
5' - GCT AGT TAT TGC TCA GCG G - 3'  
5' - ATT TAG GTG ACA CTA TAG - 3'  
5' - AAT TAA CCC TCA CTA AAG GG - 3'  
5' - CGC AAA TGG GCG GTA GGC GTG -3'  
5' - CTA CAA ACT CTT CCT GTT AGT TAG -3'  
5' - ATG GCT CAT AAC ACC CCT TG -3'  
5' - CCG GGA GCT GCA TGT GTC AGA GG -3'  
5' - GGG CTG GCA AGC CAC GTT TGG TG -3'  
5' - CTA GCA AAA TAG GCT GTC CC -3'  
5' - CTT TAT GTT TTT GGC GTC TTC CA -3'  
5' - GCA AAT GGC ATT CTG ACA TCC -3'  
5' - GAC TGG TTC CAA TTG ACA AGC -3'  
5' - CTA GAA GGC ACA GTC GAG GC -3'  
5' - ATG CGT CCG GCG TAG A -3'  
5' - GGT CCT TCT TGA GTT TGT AAC -3'  
5' - CCA TCT AAT TCA ACA AGA ATT GGG ACA AC -3'  
5' - ACC CTA ACT GAC ACA CAT TCC -3'  
5' - CTT TAT CCA GCC CTC AC -3'  
5' - GAG GTA TAT ATT AAT GTA TCG -3'  
5' - GAT TTA ATC TGT ATC AGG -3'

Contact us now at

[oligos@aitbiotech.com](mailto:oligos@aitbiotech.com)

**AIT**biotech