

01/05/2025 to 31/08/2025

Your Research, Our Priority - Discover AITbiotech Promotions Today



Offer is valid while stock lasts!



PROMO 1

All *abWare* consumables


BUY 10 GET 4 FREE

Mix and match is allowed. Free items need to be lowest in value.




Micro Tubes

	Part No.	Description
	800013	<i>abWare</i> ™ 1.5ml Microcentrifuge Tube, Conical-Bottom, Graduated, Lock Cap, Pyrogen/DNase/RNase-free, 500/Pack

PCR Consumables

	Part No.	Description
	800018	<i>abWare</i> ™ 0.2ml PCR Tube with attached Flat Cap, Clear, Polypropylene, Autoclavable, Pyrogen/DNase/RNase-free, 1000/Pack
	800019	<i>abWare</i> ™ 0.2ml 8-tube PCR Strip with attached Flat Cap, Clear, Polypropylene, Autoclavable, Pyrogen/DNase/RNase-free, 120/Pack
	800033	<i>abWare</i> ™ 0.1ml 8-tube qPCR Strip & Flat 8-cap Strip, Clear, Pyrogen/DNase/RNase-free, Autoclavable, Polypropylene, 125/Pack
	800034	<i>abWare</i> ™ 0.2ml 8-tube qPCR Strip & Flat 8-cap Strip, Clear, Pyrogen/DNase/RNase-free, Autoclavable, Polypropylene, 125/Pack
	800035	<i>abWare</i> ™ 0.1ml 8-tube qPCR Strip & Flat 8-cap Strip, White, Pyrogen/DNase/RNase-free, Autoclavable, Polypropylene, 125/Pack
	800036	<i>abWare</i> ™ 0.2ml 8-tube qPCR Strip & Flat 8-cap Strip, White, Pyrogen/DNase/RNase-free, Autoclavable, Polypropylene, 125/Pack
	800045	<i>abWare</i> ™ 0.1ml 96-well qPCR Plate, No Skirt, Clear, Pyrogen/DNase/RNase Free, Autoclavable, Polypropylene, H12-Cut, 25/Pack
	800046	<i>abWare</i> ™ 0.1ml 96-well qPCR Plate, Semi Skirt, Clear, Pyrogen/DNase/RNase Free, Autoclavable, Polypropylene, A1-Cut, 25/Pack
	800047	<i>abWare</i> ™ 0.2ml 96-well qPCR Plate, No Skirt, Clear, Pyrogen/DNase/RNase Free, Autoclavable, Polypropylene, H1-Cut, 25/Pack
	800048	<i>abWare</i> ™ 0.2ml 96-well qPCR Plate, Semi Skirt, Clear, Pyrogen/DNase/RNase Free, Autoclavable, Polypropylene, A12-Cut, 25/Pack

Microtube/Tube Racks

	Part No.	Description
	800109A	<i>abWare</i> ™ 24-Well Reversible Microtube Rack for 0.5ml or 1.5/2ml, Assorted Colours, 1/Pack
	800118	<i>abWare</i> ™ 60-Well Reversible Microtube Rack for 0.5ml or 1.5/2ml, White, 1/Pack
	800120V	<i>abWare</i> ™ 80-Well Reversible Microtube Rack for 1.5/2ml or 60 x 0.2ml, Violet, 1/Pack
	800112	<i>abWare</i> ™ 96-well Reversible Microtube Rack for 0.5ml or 1.5/2ml, Yellow, 1/Pack
	800130B 800130G 800130O	<i>abWare</i> ™ 4-Way Tube Rack, Blue, 2/Pack <i>abWare</i> ™ 4-Way Tube Rack, Green, 2/Pack <i>abWare</i> ™ 4-Way Tube Rack, Orange, 2/Pack
	800131A	<i>abWare</i> ™ Multirack for 0.2ml - 50ml Tubes, Assorted Colours, 1/Pack
	800144B 800144O 800144G	<i>abWare</i> ™ 50-Well Multi-purpose Tube Rack, Diameter: 18mm, Blue, 1/Pack <i>abWare</i> ™ 50-Well Multi-purpose Tube Rack, Diameter: 18mm, Orange, 1/Pack <i>abWare</i> ™ 50-Well Multi-purpose Tube Rack, Diameter: 18mm, Green, 1/Pack
	800150A	<i>abWare</i> ™ Snap-Together Tube Rack, Diameter: 21mm x 40 tubes (250 x 105 x 65mm), Assorted Colour, 1/Pack
	800151A	<i>abWare</i> ™ Snap-Together Tube Rack, Diameter: 16mm x 60 tubes (248 x 105 x 65mm), Assorted Colour, 1/Pack

Freezer Storage Boxes

	Part No.	Description
	800208	<i>abWare™</i> 50-well Freezer Box with Cover, Transparent, 1/Pack
	800203A	<i>abWare™</i> 100-well Freezer Storage Box, Polypropylene, Assorted colors, 1/Pack
	800248	<i>abWare™</i> 100-well Freezer Storage Box, with Gridlines on Cover, Polypropylene, White, 1/Pack

Microscope Slide Storage Box/Folder

	Part No.	Description
	800240B	<i>abWare™</i> 100-place Slide Storage Box with buckle, Blue, 1/Pack
	800254	<i>abWare™</i> 100-place Slide Storage Box with buckle, Beige, 1/Pack
	800254R	<i>abWare™</i> 100-place Slide Storage Box with buckle, Red, 1/Pack
	800254O	<i>abWare™</i> 100-place Slide Storage Box with buckle, Orange, 1/Pack
	800255	<i>abWare™</i> Cardboard Slide Folder for 20 Microscope Slides, 5/Pack

Others

	Part No.	Description
	800154	<i>abStand™</i> , linear (For <i>abCurate™</i> Pipettes), White, 1/Pack
	800451C	<i>abWare™</i> 60ml Reagent Reservoir, V-shaped bottom with lid, Polypropylene, Autoclavable, DNase/RNase Free, Pathogen-free, 5/Pack, 25pcs/Case



3 pcs of *abCurate™* Pipettes

**choose from 0.5-10ul, 10-100ul or 100-1000ul*

1 UNIT of *abStand™* (800154)



1 PACK of *abWare™*
Microcentrifuge Tubes

**Choose from 1.5ml (800013) or 2.0ml (800014)*

While stock lasts

abCurate™ Single Channel Pipettes

abCurate™ Pipettes are our newly launched powerful single channel manual air displacement pipettes with high accurate and precise aspiration and dispensing of liquid volumes. Qualified accordingly to ISO 8655/DIN 12650, these reliable pipetting tools offer the following features:

- Smooth and low operating force reduce the risk of RSI (Repetitive Strain injury).
- Presence of lockable adjustment knob to prevent accidental adjustment of set volume during pipetting.
- Clear 4-digit volume display window.
- Compatible with a wide range of pipette tips.
- Fully-Autoclavable at 121°C for 20mins.

Order Information

Part. No	Description
800902	<i>abCurate™</i> Pipette, Fully Autoclavable, Single-channel, 0.5-10 µl
800904	<i>abCurate™</i> Pipette, Fully Autoclavable, Single-channel, 10-100 µl
800907	<i>abCurate™</i> Pipette, Fully Autoclavable, Single-channel, 100-1000 µl

For more information

📍 25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477
 ☎ (65) 6778 6822 ✉ sales@aitbiotech.com 🌐 www.aitbiotech.com



PROMO 3



abSPIN™ Fast



**2nd UNIT
AT
50% OFF**

While stock lasts



abSPIN™ Fast Mini Centrifuge

abSPIN Fast Mini Centrifuge is a powerful and reliable fast spin centrifuge that is ideal for molecular biology separation, microfiltration, and quick spin of microtubes and PCR tubes.



Adjustable speed from
2000rpm to 12,000rpm.



Economical



Compact in size



Rotors are Interchangeable.
No additional tool is required
to change its rotor.


Order Information

Part. No	Description
800505	abSPIN FastMini Centrifuge complete with 8 x 1.5/2.0ml rotor and PCR Strip Tubes Rotor, 4 x 8-Tube PCR Strips, 8 x 0.5ml and 8 x 0.2ml adapters included. Do not require any tools for changing rotors.

For more information

📍 25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477

☎ (65) 6778 6822 ✉ sales@aitbiotech.com 🌐 www.aitbiotech.com



PROMO 4



abCycler™




SPECIAL PRICE



**1 PACK of abWare™
PCR Rack (800126)**

**1 PACK of abWare™ 8-tube
qPCR Strip (800034)**

While stock lasts



abCycler™ PCR Thermal Cycler

The *abCycler* is a thermal cycler intended for nucleic acid amplification through reliable and accurate Polymerase Chain Reaction (PCR). With its precise temperature controls and rapid temperature ramping, it proves to be valuable in different scientific applications, including the vast field of molecular biology, generic analysis, medical and even forensic science. With a built-in software, it operates as a standalone system thus saving space and energy at the same time.

Product Information

Sample Capacity	12 x 0.2ml 8-tube PCR strips / 96 x 0.2ml PCR tubes / 1 PCR plate 8 x 12 (No Skirt, Semi Skirt, Full Skirt)
Thermoblock	Aluminium
Temperature Range	0 – 100°C
Max. heating rate	5.5°C / sec
Max. cooling rate	4.5°C / sec
Gradient Range	1 – 40°C
Gradient block	12 columns
Gradient temperature range	35 – 100°C

Order Information

Part. No	Description
800710	abCycler PCR Thermal Cycler for 96 x 0.2ml, with Gradient

For more information

📍 25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477

☎ (65) 6778 6822 ✉ sales@aitbiotech.com 🌐 www.aitbiotech.com



TRADE-IN PROMOTION

Upgrade Your Spin – Trade In & Save Big!

 **Limited Units are Available!** 

 **Trade-in is valid for any brand of microcentrifuge.** 



abSPIN™ Microcentrifuge 24



For more information with regards to the product, please refer to the QR code.



abSPIN™ Refrigerated Microcentrifuge 24R



For more information with regards to the product, please refer to the QR code.

Order Information	
Part. No	Description
800506	abSPIN Microcentrifuge 24 complete with 1 x 24 x 1.5 / 2.0ml Fixed-angle Rotor 24 x 0.2ml PCR Tubes adapters and 24 x 0.5ml Tubes adapters
800507	abSPIN Refrigerated Microcentrifuge 24R complete with 1 x 24 x 1.5 / 2.0ml Fixed-angle Rotor 24 x 0.2ml PCR Tubes adapters and 24 x 0.5ml Tubes adapters

For more information

📍 25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477

☎ (65) 6778 6822 ✉ sales@aitbiotech.com 🌐 www.aitbiotech.com



PROMO 6



Riverstones Gloves



**BUY 10
GET 2
FREE**

While stock lasts



Riverstones Gloves



***Mix and match is allowed. Free items need to be lowest in value.**

Order Information	
Part. No	Description
RSNE80400	9 Blue Nitrile Examination Glove 100% Nitrile (Acrylonitrile-butadiene), Powder Free, Non-sterile, 4.0 mil, Blue, Finger-textured. Size: X-Small (100 pieces/box, 10 boxes/Carton)
RSNE80401	9 Blue Nitrile Examination Glove 100% Nitrile (Acrylonitrile-butadiene), Powder Free, Non-sterile, 4.0 mil, Blue, Finger-textured. Size: Small (100 pieces/box, 10 boxes/carton)
RSNE80402	9 Blue Nitrile Examination Glove 100% Nitrile (Acrylonitrile-butadiene), Powder Free, Non-sterile, 4.0 mil, Blue, Finger-textured. Size: Medium (100 pieces/box, 10 boxes/carton)
RSNE80403	9 Blue Nitrile Examination Glove 100% Nitrile (Acrylonitrile-butadiene), Powder Free, Non-sterile, 4.0 mil, Blue, Finger-textured. Size: Large (100 pieces/box, 10 boxes/carton)
RSNE80404	9 Blue Nitrile Examination Glove 100% Nitrile (Acrylonitrile-butadiene), Powder Free, Non-sterile, 4.0 mil, Blue, Finger-textured. Size: X-Large (100 pieces/box, 10 boxes/carton)

For more information

📍 25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477

☎ (65) 6778 6822 ✉ sales@aitbiotech.com 🌐 www.aitbiotech.com



PROMO 7



Low Retention
Filter Tips



**BUY 5
GET 1
FREE**

While stock lasts



Nerbe plus Low Retention Filter Tips



***Mix and match is allowed. Free items need to be lowest in value.**

Order Information	
Part. No	Description
06-613-5300	Filter tip PP, premium surface, 0.1-10ul, super slim, Gilson P2/P10 and Eppendorf type, transparent, moulded rings, np premium bio, sterile R, CE/IVD (10 boxes of 96 tips/pack)
06-642-5300	Filter tip PP, premium surface, 0-100ul, transparent, moulded rings, np premium bio, sterile R, CE/IVD (10 boxes of 96 tips/pack)
06-672-5300	Filter tip PP, premium surface, 0-300ul, natural, moulded rings, np premium bio, sterile R, CE/IVD (10 boxes of 96 tips/pack)
06-695-5300	Filter tip PP, premium surface, 100-1250ul, transparent, moulded rings, np premium bio, sterile R, CE/IVD (6 boxes of 96 tips/pack)

For more information

📍 25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477

☎ (65) 6778 6822 ✉ sales@aitbiotech.com 🌐 www.aitbiotech.com

