



Available via Catalogue Purchase
(Valid from 1/10/2024 to 31/12/2024)



3 pcs of
abCurate™ Pipettes

*choose from 0.5-10ul, 10-100ul or 100-1000ul



1 unit of
abStand™
(800154)

While stock lasts



abCurate™ Single Channel Pipettes

abCurate™ Pipettes are our newly launched powerful single channel manual air displacement pipettes with high accurate and precise aspiration and dispensing of liquid volumes. Qualified accordingly to ISO 8655/DIN 12650, these reliable pipetting tools offer the following features:

- Smooth and low operating force reduce the risk of RSI (Repetitive Strain injury).
- Presence of lockable adjustment knob to prevent accidental adjustment of set volume during pipetting.
- Clear 4-digit volume display window.
- Compatible with a wide range of pipette tips.
- Fully-Autoclavable at 121°C for 20mins.



For more information with regards to the product, please refer to the QR code.

Order Information

| Part. No | Description |
|----------|--|
| 800902 | abCurate™ Pipette, Fully Autoclavable, Single-channel, 0.5-10 µl |
| 800904 | abCurate™ Pipette, Fully Autoclavable, Single-channel, 10-100 µl |
| 800907 | abCurate™ Pipette, Fully Autoclavable, Single-channel, 100-1000 µl |

For more information

25 Pandan Crescent, TIC Tech Centre, #05-15, Singapore 128477

(65) 6778 6822 | sales@aitbiotech.com | www.aitbiotech.com

